


## CuroCell® iA Automatic & CuroCell® iA Manual

Air mattress systems

Instructions for use item number: 95-001456-EN0000



7331345178875



# Table of contents

<b>1 Introduction</b> .....	<b>3</b>
1.1 General information .....	3
1.2 Intended purpose.....	3
1.3 Intended user.....	3
1.4 Intended use environment.....	3
1.5 Indications .....	3
1.6 Contraindications.....	3
1.7 Clinical benefit .....	4
<b>2 Assembly and installation</b> .....	<b>4</b>
<b>3 Common operations</b> .....	<b>6</b>
3.1 CPR (Cardiopulmonary resuscitation) ...	6
3.2 Restart .....	6
3.3 Maximum pressure notification.....	6
3.4 Power failure.....	6
<b>4 Operation CuroCell® iA Automatic</b> .....	<b>7</b>
4.1 Automatic setting .....	9
4.2 Programs .....	9
4.3 Maximum pressure (caring mode).....	9
4.4 Panel lock.....	9
4.5 Comfort settings .....	9
4.6 Hand check (function control).....	10
4.7 CPR (Cardiopulmonary resuscitation) .	10
4.8 Restart .....	10
4.9 Power failure.....	10
4.10 Notifications .....	10
4.11 Maximum pressure notifications.....	10
4.12 Table of notifications.....	10
<b>5 Operation CuroCell® iA Manual</b> .....	<b>12</b>
5.1 Manual setting.....	14
5.2 Programs .....	14
5.3 Maximum pressure (caring mode).....	14
5.4 Seating function.....	14
5.5 Sitting positioning in bed .....	14
5.6 Panel lock.....	14

5.7 Hand check (function control) .....	15
5.8 CPR (Cardiopulmonary resuscitation) ..	15
5.9 Restart.....	15
5.10 Power failure .....	15
5.11 Notifications .....	15
5.12 Maximum pressure notification .....	15
5.13 Table of notifications.....	16
<b>6 Product description</b> .....	<b>17</b>
6.1 Product combinations.....	17
6.2 Control unit (iA Automatic/iA Manual)17	
6.3 Mattresses.....	19
<b>7 Reuse and cleaning</b> .....	<b>22</b>
7.1 Cleaning and disinfection .....	22
7.2 Reconditioning.....	22
<b>8 Storage</b> .....	<b>23</b>
<b>9 Service and maintenance</b> .....	<b>24</b>
9.1 Service and maintenance chart .....	24
9.2 Replacing the air filter .....	25
<b>10 Troubleshooting</b> .....	<b>25</b>
<b>11 Technical specification</b> .....	<b>26</b>
11.1 Standards .....	27
11.2 Symbol key .....	28
<b>12 Other information</b> .....	<b>29</b>
12.1 Recommended lifetime.....	29
12.2 Disassembly and recycling.....	29
12.3 Returns .....	29



## WARNINGS AND SAFETY PRECAUTIONS

**Read all instructions before use or repair.**

**WARNING! To minimize the risk of fire, personal injury and equipment/property damage adhere to the following instructions:**

1. The product must only be installed and used for its intended purpose according to the instructions in this manual and/or other documentation from Care of Sweden. The product must not be combined, assembled or repaired with parts (e.g. control unit and mattress), accessories or spare parts other than those described in this manual or other documentation from Care of Sweden. The product must not be modified in any way.
2. The product must be placed and used so that it does not become trapped or damaged. Be particularly aware of trapping damage when using side rails. Always make sure that the mattress is the correct size for the bed.
3. Regularly check product functionality by performing a hand check.
4. When the product is used for patients needing special supervision, such as children, continuous monitoring is required.
5. The mattress is protected by a hygiene cover; avoid using multiple hygiene covers as this can affect the vapor permeability of the mattress.
6. The hygiene cover does not allow liquid or air to penetrate but is vapour permeable. Make sure that the patient is positioned correctly to avoid the risk of suffocation.
7. Be careful with sharp objects to prevent damage to the hygiene cover.
8. Do not open the control unit housing – risk for electric shock. Servicing and maintenance must be performed by Care of Sweden or one of its authorized service technicians.
9. Route the power cable to the control unit carefully to avoid tripping. Also make sure that the patient is lying correctly on the mattress according to the instructions and use a cable holder if possible.
10. To avoid the risk of strangulation, make sure that the cable and tubes are routed to prevent someone getting caught up in them.
11. Do not use the product in bathrooms or other areas where there is a risk of the control unit coming into contact with water or other liquids. Except for specified cleaning, never handle a product that has come into contact with water/liquid. Pull the plug out of the electrical socket immediately and send the product to an authorized service technician for servicing.
12. Do not use close to or in contact with fire sources/hot surfaces, such as fire, burning cigarettes, hot lamps, heating fans or heating stoves/open fires as this could damage the product.
13. Do not store or use the product in direct sunlight. The product may be damaged by the elevated temperature and UV light.

14. Strong magnetic fields or wireless communication equipment (e.g., wireless home network products, mobile phones, walkie-talkies, cordless phones and their base stations, radio transmitters, etc.) may affect the product's functionality and should be kept at a distance of at least 1 meter from the control unit.
15. Never use the product if the power cable, plug of the control unit or power supply housing is defective, if the control unit housing is damaged, or if it is not functioning properly. Contact an authorized service technician for examination and repair.
16. Never connect anything other than the Care of Sweden supplied power supply to the control unit power cable connector.
17. If the hygiene cover is equipped with side handles, these are intended for managing or relocating the mattress. Do not use the handles to lift the mattress with the patient lying on it. All other use takes place under your own liability and is not covered by the product warranty.
18. To prevent the power supply from being pulled out, exercise caution when there are children and pets in the environment around the equipment.
19. Use of this product adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this product and the other equipment should be observed to verify that they are operating normally.
20. To minimize the risk of wounds occurring on the feet, make sure that the patient doesn't come into contact with the hangers of the control unit.


# 1 Introduction

These air mattress systems may be used as an aid to prevent and treat pressure ulcers/pressure injuries.

**CuroCell® iA Automatic** is a control unit that sets the air pressure in the mattress based on the patient's length, weight, and position without need for manual settings. It offers the operator the option to adjust program and comfort.

**CuroCell® iA Manual** is a control unit where the operator must manually adjust the pressure in the mattress based on the weight of the patient.

**CuroCell® iA Automatic** and **CuroCell® iA Manual** are compatible with six different mattresses: CuroCell® Ci10, CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci10 PRO, CuroCell® Ci17 PRO and CuroCell® Ci20 PRO. See more information about the mattresses in section 6.3.

 Always read the instructions for use prior to use.

## 1.1 General information

The system is a medical device with CE marking in accordance with MDR (EU) 2017/745. According to this regulation the manufacturer is required to report all accidents or incidents involving the products. All information involving accidents or incidents relating to our products, shall be reported immediately to Care of Sweden.

## 1.2 Intended purpose

The mattress system consists of a control unit and a mattress and is intended to be used for prevention and as an aid in the treatment of pressure ulcers/pressure injuries (PU/PI).

## 1.3 Intended user

The mattress system is intended to be used by all kind of patients, including lay persons. Note that the patient and operator could be the same person. The mattresses are intended for use by patients of a recommended minimum length of 120 cm.

The specifications for weight are listed in the following table.

Mattress	Recommended patient weight
CuroCell® Ci10	≤ 170 kg
CuroCell® Ci17	≤ 200 kg
CuroCell® Ci20	≤ 230 kg
CuroCell® Ci10 PRO	≤ 180 kg
CuroCell® Ci17 PRO	≤ 210 kg
CuroCell® Ci20 PRO	≤ 240 kg

## 1.4 Intended use environment

The mattress system is intended to be used in all kinds of health care environments, including home care.

## 1.5 Indications

Suitable for a wide range of patients with increased risk for pressure ulcers/injuries, including those with superficial ulcers, up to category IV in association with an individualized plan of care.

## 1.6 Contraindications

There are no known contraindications. It is necessary for the prescriber of the mattress to make an individual assessment of the patient and decide whether the characteristics of the mattress is beneficial for the patient based on diagnosis.

## 1.7 Clinical benefit

The clinical benefits for CuroCell® iA Automatic and CuroCell® iA Manual together with any of the mattresses included in this instruction for use are:

- Prevention and treatment of pressure ulcers/pressure injuries up to and including category IV.
- Reduction of shear forces.
- A safe, comfortable pressure redistribution mattress system which is easy to handle.
- Silent running control units.

### Note!


For certain patients, e.g., amputees, the recommended length measurement may not be reached. Patients in these groups may require other settings as the entire surface is not under load. For function controls, see section 4 or 5 depending on which control unit you have. The mattress may be inappropriate for use during x-ray examinations because of the risk of blurred images or artefacts that may lead to diagnostic errors. In the event of a power loss or similar, the mattress will retain air for at least 12 hours.

## 2 Assembly and installation

Check that no parts are damaged. If damage is found, contact Care of Sweden or your local distributor before using the product. Do not use sharp objects when unpacking as it might damage the product.

**When using an overlay mattress (CuroCell® Ci10 or CuroCell® Ci10 PRO):**


**1a.** Place the overlay mattress on the base mattress. Secure the mattress to the base using the 4 straps on the corners of the mattress.



**When using a full replacement mattress (CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci17 PRO or CuroCell® Ci20 PRO):**

**1b.** Place the mattress on the bed base.

**1c.** When using CuroCell® Ci17 PRO or CuroCell® Ci20 PRO: Secure the mattress to the bed using the



fastening straps on the underside of the mattress. If the mattress is used on an adjustable bed, the 4 straps at the head end shall be fastened at the moving part of the bed's head end. The 2 straps at the foot end of the bed shall be fastened at the moving part of the foot end.


**1d.** When using CuroCell® Ci17 or CuroCell® Ci20 the mattress is kept in position by the Grip-lock cover.

### Note!

Make sure that the mattress is the correct size for the bed. Check the cells and press studs to ensure they are correctly assembled. Make the bed with sheets. The mattress should be used lying in the lengthwise direction on the mattress, with feet at the end, marked with the feet symbol.




**2.** Hang the control unit on the foot end of the bed or place it on a level, steady surface. Check that the switch on the side of the control unit is set to '0' (off).




**3.** Place the power cable so there is no risk of stumbling over it, running over it with the bed wheels, or getting it jammed when raising or lowering the bed. Put the power cable into the control unit.


4. Open the lid on the air tube connector (marked CPR) and connect it to the side of the control unit.



5. A click is heard and felt when correctly connected. Secure that both sides of the connection are closed.



6. Plug the power supply into an approved and easily accessible electrical socket (100– 240 V).



7. Check that the power cable has been correctly connected to the control unit and that the correct power supply has been used. (See section 11 Technical specification). The correct P/N is shown on the label on the power supply. The power supply is part of the equipment and may not be replaced.

If the control unit has been stored in its minimum or maximum storage temperature (-25°C to +70°C), wait at least 1 hour before starting it. This time is based on an ambient temperature of 20°C.

**Note!**  
Do not hold the 12V plug on the power supply while touching the patient.

### 3 Common operations


Following operations apply to CuroCell® iA Automatic and CuroCell® iA Manual regardless of which mattress that is used.

#### 3.1 CPR (Cardiopulmonary resuscitation)

In case of an emergency where CPR (Cardiopulmonary resuscitation) is necessary, remove the connection (marked 'CPR') from the control unit and leave the lid open to empty the mattress of air quickly.

#### 3.2 Restart

If a restart is required, set the On/Off switch on the side of the control unit to 0 (Off). Wait for approx. 10 seconds and restart the control unit.



#### 3.3 Maximum pressure notification



When the function Maximum pressure has been used for a long time, the Maximum pressure diode will blink.

If the use is intentional, ignore the notification.

#### 3.4 Power failure

In the event of a power failure, the valves will open automatically and even out the air pressure in the mattress. The mattress will retain air for at least 12 hours. Perform a hand check to make sure the pressure of the mattress is not too hard or too soft. See 4.6 for CuroCell® iA Automatic and 5.7 for CuroCell® iA Manual for more information about hand check.



## 4 Operation CuroCell® iA Automatic


Following operations apply only to CuroCell® iA Automatic regardless of which mattress that is used. Read the label of the control unit carefully to make sure that you know which control unit you have.

Operation CuroCell® iA Automatic (how to get started)
1. Set the On/Off switch on the side of the control unit to position 1 (On). See figure 1.
2. The mattress starts to inflate. This takes about 20-40 minutes depending on the size of the mattress. While the mattress is inflating, the "mute the information signal" and "information signal" diodes light up in orange. When these diodes goes out, the patient can be placed on the mattress.
3. The Gentle Alternating Low Pressure (GALP) setting is pre-set. If a different setting is desired, it can be selected when the patient is placed on the mattress.
4. The control unit sets the inner pressure of the mattress according to the weight, length, and position of the patient. This takes about 20-30 minutes. During this time, the diode above the selected program flashes. When the diode stops flashing, the mattress is set and the inner pressure has adjusted according to the patient.
5. Perform a hand check to ensure that the settings are correct.

### Note!










Each time the system starts up, it will operate as follows:

- When using a mattress without an air safety mattress (CuroCell® Ci10 or CuroCell® Ci10 PRO), the control unit will completely inflate the main cells and then perform an automatic setting.
- When using a mattress with an air safety mattress (CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci17 PRO or CuroCell® Ci20 PRO) the control unit always begins by completely inflating the safety mattress, followed by the other cells, and then perform an automatic setting.
- The mattress must be inflated before the patient can lie on the mattress.
- During the automatic setting, try to avoid larger movements on the mattress. Otherwise, you will get a notification that the desired value could not be reached within the time limit and the weighting must start over.
- Once the automatic setting is complete, the control unit switches to a basic setting of Gentle Alternating Low Pressure (GALP) (when used for the first time) or to the previous setting. For more information see section 4.1.



**Figure 1.** On/Off switch.



Button	Function
	Mute the information signal
	Panel lock
	Gentle Alternating Low Pressure (GALP)
	Constant Low Pressure (CLP)
	Maximum pressure (caring mode)
	Information signal
	Incorrect connection of the air connector (CPR)
	Cycle time settings (10, 15, 20, 25 minutes). The diodes are also used for error notifications.
	Comfort settings

## 4.1 Automatic setting

The mattress system independently and without manual adjustment controls the inner pressure of the mattress according to the weight, height, and position of the patient. No manual action needs to be performed to adjust the inner pressure of the mattress. This function works as follows:

1. At start-up, automatic setting of the inner pressure of the mattress is always carried out according to the weight and height of the patient.
2. If the patient moves noticeably or changes their position, the system will independently control the inner pressure of the mattress.
3. The system performs an automatic setting at fixed intervals even if no significant changes have occurred.

After automatic setting of the mattress inner pressure, the system returns to the previously selected program. At start-up, Gentle Alternating Low Pressure is always pre-set.

## 4.2 Programs

There are two programs to choose from:

**Constant Low Pressure (CLP)** means that the air pressure in all the air cells is the same in the whole mattress.

**Gentle Alternating Low Pressure (GALP)** means that the air pressures in the air cells are different and alternates regularly after chosen cycle periods.

Choose program by pushing the button for the program. We recommend Gentle Alternating Low Pressure (GALP) which is also the preset mode.



1. **Constant Low Pressure (CLP).**  
No cycle period is needed.



2. **Gentle Alternating Low Pressure (GALP).** The cycle period can be changed according to the patient's needs and requirements. Choose between 10, 15, 20 or 25 minutes. The longer the cycle period, the slower the alternations. We recommend a basic setting of 10 minutes.

## 4.3 Maximum pressure (caring mode)



With this function, the entire mattress is inflated and provides firmed support. This function reverts to the previous setting after 20 minutes. The function should be used when caring the patient, shifting the patient's position, or moving the patient in or out of bed.

## 4.4 Panel lock



Press the Panel lock button to lock or unlock the control panel. The button indicates when the panel has been locked. The screen locks automatically if left untouched for five minutes. This is to prevent the settings being changed accidentally. To unlock, press the button for 2 seconds.

## 4.5 Comfort settings



The pressure can be increased in two steps depending on the patient's comfort requirements. This increase is made based on the automatic setting in 4.1. The selected setting is shown by a green light.

### Note!

- When only parts of the mattress are under load, for example, in the case of amputees – it may be necessary to raise the setting using the comfort settings.

#### 4.6 Hand check (function control)

Hand check is performed to ensure that the mattress system works properly. Hand check should be performed regularly; for CuroCell® iA Automatic we recommend once every eight hours as well as after installation of the system.

##### Note!

- Make sure that the mattress system is filled, which is shown by a green light from the diode, before performing a hand check.
- How to perform a hand check depends on which mattress is used - read the label of the mattress carefully to know which mattress you have.

#### When using an overlay mattress system

(CuroCell® Ci10 or CuroCell® Ci10 PRO)

1a. Open the cover and identify a cell with less air at the patient's sacrum. Insert a hand, with the palm facing up, between the overlay mattress and the underlying mattress. The hand is inserted beneath the patient's sacrum (center of mattress).

#### When using a full replacement mattress

(CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci17 PRO or CuroCell® Ci20 PRO)

1b. Open the cover and identify a cell with less air at the patient's sacrum. Insert a hand, with the palm facing up, between the top cells and the underlying safety mattress. The hand is inserted beneath the patient's sacrum (center of mattress).

2. Ensure there is a gap between the patient and the underlying mattress so that the patient does not 'bottom out'.

3. If you can feel the patient's sacrum resting in the palm of your hand, the gap is too small. See section 10 'Troubleshooting'.

#### 4.7 CPR (Cardiopulmonary resuscitation)

See 3.1

#### 4.8 Restart

See 3.2


#### 4.9 Power failure

See 3.4

#### 4.10 Notifications



Different notifications exist based on how serious the warning is. With a malfunction or an error, a notification will be given by a flashing warning triangle. To mute the warning signal, press the mute button.



When a notification occurs, the current cycle time diode will turn off and a notification code is shown on the four different cycle time diodes (10, 15, 20, 25). To read the cycle time during the error notification, unlock the control panel.














#### 4.11 Maximum pressure notifications

See 3.3

#### 4.12 Table of notifications

Information about each notification is shown in the notification table on the next page.

- Some of the notifications are both audible and visual. The error code will be displayed until the error has been rectified. If the mute button is pressed the audible warning will cease for a period of 5 minutes and will return until the error has been rectified.
- Some of the notifications have no audible alarm. The error code is shown until the system is restarted.

Notification (audible and visual)	Description and troubleshooting
	High temperature. Valves and compressors are turned off. If the control unit is in direct sunlight, relocate it. If problem persists, contact an authorized service technician.
	Default settings are not completed. Contact an authorized service technician.
	Incorrect input voltage. Make sure the correct power supply is used, otherwise contact an authorized service technician.
	The display is the probable cause. Contact an authorized service technician.
	Low pressure. Check the CPR, mattress, air tubes and air filter. If problem persists, contact an authorized service technician.
	High pressure. The pressure cannot be reduced to the desired value within the time limit. If problem persists, contact an authorized service technician.
	The automatic setting has been restarted too many times during the automatic setting period. Alarm can be caused by too much movement. Restart the control unit, have the user lay still during start-up. If problem persists, contact an authorized service technician.
	Mattress error read. The system cannot detect a mattress connected to the control unit. Check the CPR. If problem persists, contact an authorized service technician.
	The pressure is not increasing fast enough during operation. Check the CPR, air tubes and air filter. If problem persists, contact an authorized service technician.
Notification (only visual)	Description and troubleshooting
	Leakage in cell one, cell two or valve. Check the CPR, mattress, air tubes and air filter. If problem persists, contact an authorized service technician.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Contact an authorized service technician. Information to service technicians: this notification shows a leakage in the blue section.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Contact an authorized service technician. Information to service technicians: this notification shows a leakage in the green section.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Contact an authorized service technician. Information to service technicians: this notification shows a leakage in the red section.

## 5 Operation CuroCell® iA Manual


Following operations apply only to CuroCell® iA Manual regardless of which mattress that is used. Read the label of the control unit carefully to make sure that you know which control unit you have.

Operation CuroCell® iA Manual (how to get started)
1. Set the On/Off switch on the side of the control unit to position 1 (On). See figure 1.
2. Use the weight settings to set a suitable air pressure.
3. The mattress starts to inflate. This takes about 20-40 minutes depending on the size of the mattress. While the mattress is inflating, the “mute the information signal” and “information signal” diode light up in orange. When these diodes go out, the patient can be placed on the mattress.
4. The Gentle Alternating Low Pressure (GALP) setting is pre-set. If a different setting is desired, it can be selected when the patient is placed on the mattress.
5. When the patient is placed on the mattress and the desired program is selected, the inner pressure of the mattress is adjusted according to the selected weight setting. This takes about 20-30 minutes. During this time, the diode above the selected program flashes. When the diode stops flashing, the mattress is set and the inner pressure has reached the selected value.
6. Perform a hand check to ensure that the settings are correct.
7. If the mattress is too soft or too hard, adjust the weight setting.
8. Perform another hand check.

### Note!

Each time the system starts up, it will operate as follows:

- When using a mattress without an air safety mattress (CuroCell® Ci10 or CuroCell® Ci10 PRO), the control unit will adjust to the set weight.
- When using a mattress with an air safety mattress (CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci17 PRO or CuroCell® Ci20 PRO), the control unit always begins by completely inflating the safety mattress, followed by the other cells, and then start on the default mode.
- The mattress must be inflated before the patient can lie on the mattress.



**Figure 1.** On/Off switch.



Button	Function
	Mute the information signal
	Panel lock
	Gentle Alternating Low Pressure (GALP)
	Constant Low Pressure (CLP)
	Maximum pressure (caring mode)
	Patient weight settings
	Information signal
	Incorrect connection of the air connector (CPR)
	Seating function
	Cycle time settings (10, 15, 20, 25 minutes). The diodes are also used for error notifications.

## 5.1 Manual setting

At start-up, the inner pressure of the mattress must be set manually based on the weight and height of the patient. The mattress system maintains the pre-set inner pressure regardless of movement and position changes. This means that when the patient changes position, for example, the mattress's inner pressure must be adjusted manually. The weight settings on the control unit are used to change the inner pressure.

## 5.2 Programs

There are two programs to choose from:

**Constant Low Pressure (CLP)** means that the air pressure in all of the air cells is the same in the whole mattress.

### Gentle Alternating Low Pressure (GALP)

means that the air pressures in the air cells are different and alternates regularly after chosen cycle periods.

Choose program by pushing the button for the program. We recommend the Gentle Alternating Low Pressure mode which is also the preset mode.



### 1. Constant Low Pressure (CLP).

No cycle period is needed.



### 2. Gentle Alternating Low Pressure (GALP).

The cycle period can be

changed according to patient needs and requirements. Choose between 10, 15, 20 or 25 minutes. The longer the cycle period, the slower the alternations. We recommend a basic setting of 10 minutes.

## 5.3 Maximum pressure (caring mode)

With this function, the entire mattress is inflated and provides firmed support. This function reverts to the previous setting after approximately 20 minutes. The function should be used when caring the patient, shifting the patient's position, or moving the patient in or out of bed. When the function has been used for a long time, the diode will blink. If the use is intentional, ignore the notification.



## 5.4 Seating function



A visible and audible notification is generated when the seating function has been active for two consecutive hours, notifying that repositioning of the individual may be necessary. The control unit will also notify if an attempt is made to switch on the system while seating function is active.

## 5.5 Sitting positioning in bed

When raising the head end of the bed into a sitting position, always secure the patient's position. To ensure the product functionality, we always recommend performing a hand check, see 5.7. This function is recommended to use for short periods only. For additional support, positioning pillows can be used.

### Note!

- When using a lift to place the patient in the bed and the head end of the bed is raised, make sure that the patient is not placed too high on the mattress. Otherwise, there is a risk of shear.
- When the Gentle Alternating Low Pressure (GALP) program is used and the head end of the bed is raised, make sure that the patient and/or the mattress is not moving downwards due to the movement in the mattress. Also raise the foot end of the bed.

## 5.6 Panel lock



Press the Panel lock button to lock or unlock the control panel. The button indicates when the panel has been locked. The screen locks automatically if left untouched for five minutes. This is to prevent the settings being changed accidentally. To unlock, press the button for 2 seconds.



## 5.7 Hand check (function control)

Hand check is performed to ensure that the mattress system works properly and to ensure that the weight setting is correct. Hand check should be performed regularly; for CuroCell® iA Manual, we recommend once every eight hours as well as after installation of the system.

For manual systems hand check should also be performed when moving from lying to seating position in bed, in case of positioning changes and when changing the comfort setting (and/or the weight setting).

### Note!

- Make sure that the mattress is filled, and that the system has adapted to the patient before performing a hand check.
- How to perform a hand check depends on which mattress is used – read the label of the mattress carefully to know which mattress you have.

**When using an overlay mattress system**  
(CuroCell® Ci10 or CuroCell® Ci10 PRO):

**1a.** Open the cover and identify a cell with less air at the patient's sacrum. Insert a hand, with the palm facing up, between the overlay mattress and the underlying mattress. The hand is inserted beneath the patient's sacrum (center of mattress).

**When using a full replacement mattress**  
(CuroCell® Ci17, CuroCell® Ci20, CuroCell® Ci17 PRO or CuroCell® Ci20 PRO)

**1b.** Open the cover and identify a cell with less air at the patient's sacrum. Insert a hand, with the palm facing up, between the top cells and the underlying safety mattress. The hand is inserted beneath the patient's sacrum (center of mattress).

**2.** Ensure there is a gap between the patient and the underlying mattress so that the patient does not 'bottom out'.

**3.** If you can feel the patient's sacrum resting in the palm of your hand, the gap is too small. See section 10 'Troubleshooting'.

## 5.8 CPR (Cardiopulmonary resuscitation)

See 3.1

## 5.9 Restart

See 3.2

## 5.10 Power failure

See 3.4

## 5.11 Notifications



Different notifications exist based on how serious the warning is. With a malfunction or an error, a notification will be given by a flashing warning triangle. To mute the warning signal, press the mute button.



When a notification occurs, the current cycle time diode will turn off and a notification code is shown on the four different cycle time diodes (10, 15, 20, 25). To read the cycle time during the error notification, unlock the control panel.













## 5.12 Maximum pressure notification

See 3.3

### 5.13 Table of notifications

Information about each notification is shown in the notification table:

- Some of the notifications are both audible and visual. The error code will be displayed until the error has been rectified. If the mute button is pressed the audible warning will cease for a period of 5 minutes and will return until the error has been rectified.
- Some of the notifications have no audible alarm. The error code is shown until the system is restarted.

Notification (audible and visual)	Description and troubleshooting
	High temperature. Valves and compressors are turned off. If the control unit is in direct sunlight, relocate it. If problem persists, contact an authorized technician.
	Default settings are not completed. Contact an authorized technician.
	Incorrect input voltage. Make sure the correct power supply is used, otherwise contact an authorized technician.
	The display is the probable cause. Contact an authorized technician.
	Low pressure. Check the CPR, mattress, air tubes and air filter. If problem persists, contact an authorized technician.
	High pressure. The pressure cannot be reduced to the desired value within the time limit. If problem persists, contact an authorized technician.
	Mattress error read. The system cannot detect a mattress connected to the control unit. Check the CPR. If problem persists, contact an authorized technician.
	The pressure is not increasing fast enough during operation. Check the CPR, air tubes and air filter. If problem persists, contact an authorized technician.
Notification (only visual)	Description and troubleshooting
	Leakage in cell one, cell two or valve. Check the CPR, mattress, air tubes and air filter. If problem persists, contact an authorized technician.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Information to service technicians: this notification shows a leakage in the blue section. Contact an authorized technician.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Information to service technicians: this notification shows a leakage in the green section. Contact an authorized technician.
	Leakage in one of the sections. Secure the CPR, mattress, and connection tubes. Information to service technicians: this notification shows a leakage in the red section. Contact an authorized technician.

## 6 Product description

Make sure that you read about the right product by looking at the label of the control unit as well as the mattress.

### 6.1 Product combinations

Below is a table showing which products that are possible to combine with each other.

Control unit	Mattress	Mattress cover
CuroCell® iA Automatic CuroCell® iA Manual	CuroCell® Ci10	Cover Olivia
	CuroCell® Ci10 PRO	Cover Stone
	CuroCell® Ci17	Cover Olivia Grip-lock
	CuroCell® Ci20	Cover Stone Grip-lock
	CuroCell® Ci17 PRO	Cover Top part Olivia
	CuroCell® Ci20 PRO	Cover Top part Stone Cover Bottom part CuroCell






### 6.2 Control unit

**CuroCell® iA Automatic** is an automatic air mattress system used as an aid to prevent and treat pressure ulcers/pressure injuries. The automation means that the control unit's built-in sensors use software to adjust the inner pressure of the mattress according to the patient's weight, height, position, and change in position. This means that no manual action needs to be performed to adjust the inner pressure of the mattress to conform to the patient.

**CuroCell® iA Manual** is a manual air mattress system used as an aid to prevent and treat pressure ulcers/pressure injuries. The inner pressure of the mattress must be set manually based on the weight of the patient. Perform a hand check to make sure the pressure is correct, see section 5.7. The mattress system maintains the pre-set inner pressure regardless of movement and position changes. This means that when the patient changes position, for example, the mattress's inner pressure must be adjusted manually.

Any of these two control units are compatible with any of the mattresses in 6.3.

- 1. Control panel
- 2. Tube/CPR connection
- 3. Power switch, On/Off
- 4. Connection power cable
- 5. Air filter
- 6. Hangers
- 7. Power supply



### 6.3 Mattresses

#### CuroCell® Ci10

1. Mattress
2. Hygiene cover
3. Main cells
4. Heel cells
5. Tubing set
6. CPR (quick deflation)




#### CuroCell® Ci10 PRO

1. Mattress
2. Hygiene cover
3. Main cells
4. Cell holder (integrated in inner cover)
5. Heel cells
6. Press studs for inner cover
7. Tubing set
8. CPR (quick deflation)




### CuroCell® Ci17

- 1. Mattress
- 2. Hygiene cover - top part
- 3. Main cells
- 4. Heel cells
- 5. Tubing set
- 6. CPR (quick deflation)
- 7. Safety cells
- 8. Hygiene cover – bottom part



### CuroCell® Ci17 PRO

- 1. Mattress
- 2. Hygiene cover – top part
- 3. Main cells
- 4. Heel cells
- 5. Press studs for inner cover
- 6. Cell holder (integrated in inner cover)
- 7. Safety cells
- 8. CPR (quick deflation)
- 9. Tubing set
- 10. Hygiene cover - bottom part




### CuroCell® Ci20

1. Mattress
2. Hygiene cover – top part
3. Main cells
4. Heel cells
5. Tubing set
6. CPR (quick deflation)
7. Safety cells
8. Hygiene cover – bottom part



### CuroCell® Ci20 PRO

1. Mattress
2. Hygiene cover – top part
3. Main cells
4. Heel cells
5. Press studs for inner cover
6. Cell holder (integrated in inner cover)
7. CPR (quick deflation)
8. Tubing set
9. Safety mattress
10. Hygiene cover – bottom part



## 7 Reuse and cleaning

The product is reusable. Before reusing, it is important to follow the instructions below for cleaning, disinfection, and reconditioning. Disinfection is recommended between patients according to the instructions below.

Always follow local instructions and the instructions for use of the cleaning and disinfecting agent. Consult your hygiene manager or Care of Sweden for help and instructions in case of uncertainty.

### 7.1 Cleaning and disinfection

#### CONTROL UNIT



Wipe the control unit with a damp cloth and mild detergent. Primarily use solvent-free soap with a neutral pH value. If necessary, a disinfectant and/or cleaning agent can be used such as: alcohol with or without surfactants or oxidizing solutions such as: chlorine and/or hydrogen peroxide, concentration 1000 ppm/0,1%. In exceptional cases, a maximum concentration of 10,000 ppm/1% can be used.

If another agent is used, choose one that does not harm the exterior of the control unit.

#### INNER COVER AND MATTRESS COVER

##### Wipe off



Primarily use solvent-free soap with a neutral pH value. For a daily basis cleaning, a disinfectant and/or cleaning agent can be used such as: alcohol with/without surfactants or oxidizing solutions such as: chlorine and/or hydrogen peroxide, concentration 1000ppm/0,1%. In exceptional cases, a concentration of a maximum of 10,000ppm/1% can be used, then consider that high concentrations can shorten the life of the coating.

##### Mechanical cleaning



Covers consisting of several parts must be separated before washing.

### 7.2 Reconditioning

#### CONTROL UNIT

Clean the control unit according to section 7.1 Cleaning and disinfection – Control unit.

#### MATTRESS

Disconnect the tube connector from the control unit and remove the air from the mattress.

##### Cleaning of mattress

1. Clean all external surfaces of the mattress according to section 7.1 Cleaning and disinfection - Inner cover and Mattress cover - Wipe off. Ensure that all areas are free of dirt residues.
2. If the mattress is heavily soiled, it is recommended that the mattress is taken apart and cleaned, follow the instructions below according to points 3-5.
3. Remove the covers.
4. Wipe off the cells, tubing and the CPR module with a cleaning agent according to local instructions and the instructions for use of the cleaning and disinfecting agent.
5. When all parts are dry, assemble the mattress. If cells have become detached from the tubes, these must be put back according to drawing in section 6.3.

##### Note!

- Check the hygiene cover, cells and hoses at each cleaning. If the hygiene cover is damaged, it must be replaced. Also check the control unit, tube connections and power cord when cleaning. Damaged parts must be replaced or repaired.



### **Disinfection of mattress**

1. Disinfect all external surfaces of the mattress with disinfectant according to section 7.1 Cleaning and disinfection - Inner cover and Mattress cover - Wipe off. Ensure that all areas are free of dirt residue.
2. Allow the disinfectant to work according to the instructions from the manufacturer of the agent.
3. Let the cover dry.
4. If the mattress is heavily soiled, it is recommended that the mattress is taken apart and disinfected, follow the instructions below according to points 5-8.
5. Remove the covers.
6. Wipe the cells, tubes and the CPR module with a disinfectant.
7. Allow the disinfectant to work according to the instructions from the agent's manufacturer.
8. When all parts are dry, assemble the mattress. If cells have become detached from the tubes, they need to be put back according to the drawing in 6.3.

## **8 Storage**

It is advisable to store the mattress and control unit in the product bag (optional), original package or equivalent for protective storage. Handle the packaged product with caution. Do not place any heavy objects on top of it. For additional information about storage temperature, see section 11.

## 9 Service and maintenance

The control unit should be regularly serviced and inspected to maintain functionality and performance. Only use spare parts approved by Care of Sweden.

### 9.1 Service and maintenance chart

Service in the column to the right must only be performed by an authorized service technician.

	<b>Before every use (or every 2 weeks if used continuously)</b>	<b>After each use (between users)</b>	<b>After 12 months</b>	<b>After 5 years</b>
<b>Control unit</b>				
Visual inspection of exterior	X	X		
Cleaning of exterior		X		
Visual inspection of power cable/supply	X			
Test function on control panel		X		
Visual inspection of the connecting tubes (marked CPR), they shall be positioned correctly and not leaking		X		
Control air filter and replace if needed			X	
<b>Mattress and cover</b>				
Visual inspection of cover, no damages		X		
Test function of zipper		X		
Visual inspection of connecting tubes	X			
<b>Other</b>				
Function test (start-up)		X		
Full service				X
Change valve module				X

## 9.2 Replacing the air filter

Before any maintenance is done, make sure that the control unit is turned off. Services shall not be done while using the product.

To replace the air filter:

1. Loosen the small protective plate on the rear of the control unit using a size T10 Torx screwdriver.



2. Remove the filter from the holder.
3. Place the new filter in the holder with the pink side facing outwards. Put the protective plate back in place and secure using the screws.

If the control unit is used in a dirty environment the filter should be checked regularly.

## 10 Troubleshooting

If the problem keeps occurring, please contact Care of Sweden or your local distributor.

Problem	Solution
The control unit does not start	Check that the power supply has been connected to the mains supply. Check that the LED on the power supply is showing green.
The patient is 'bottoming out'	CuroCell® iA Automatic: Restart the control unit. See section 3.2. The control unit will initiate an automatic setting. Wait until the automatic setting is complete. Perform a hand check (see section 4.6).  CuroCell® iA Manual: Restart the control unit. See section 3.2. Perform a hand check (see section 5.7). If the gap is still too small, raise the comfort setting in stages, be careful not to raise the setting too high.
The mattress moves around	Check that the mattress is fastened to the bed frame with the straps underneath (two at the head end and two on each of the long sides).
Some cells have less air	This is normal for Gentle Alternating Low Pressure (GALP), as the air supply switches between alternating cells for a predetermined cycle period (one cycle = 10–25 minutes).
The control unit makes a noise; vibrations can be felt	Check how the control unit is hanging on the bed. Resonance can occur, in parts of the bed. Remove the control unit and listen to find out if this vibrations makes a difference. The problem may be resolved by putting the control unit on a flat, steady surface or by placing a towel between the control unit and bed.

## 11 Technical specification

CONTROL UNIT SPECIFICATION		
Model		CuroCell® iA Automatic, CuroCell® iA Manual
Input voltage		100-240 V / 50-60 Hz / 0,6 A
Output voltage		12 V DC
Power supply	Ungrounded AC outlet, electrical safety class II	Use only power supply with P/N WR9QE1500LRPCIMG3138
Power consumption		Max 18 W
Electrical classification		Class II, Type BF
Fuse		No Fuse
Mode of operation	CuroCell® iA Automatic CuroCell® iA Manual	Gentle Alternating Low Pressure, and Constant Low Pressure (CLP)
Cycle time	Gentle Alternating Low Pressure	10 min, 15 min, 20 min, 25 min
Patient pressure settings	CuroCell® iA Automatic	Automatic adjustment of patient pressure (internal air pressure) in the mattress
	CuroCell® iA Manual	Operator sets the patient pressure (internal air pressure) in the mattress according to patients' weight. Correct settings to be controlled by hand check
Dimensions (L x W x H)		11 cm x 27 cm x 15,5 cm
Weight		1,7 kg
Sound pressure level according EN ISO 11201	A-weighted emission sound pressure level LpA,eq (dB)	28 dBA (at operator position) 25 dBA (at head end) when placing the control unit on the foot end.
Sound power levels according to EN ISO 3746	A-weighted sound power level L WA (dB)	<40 dBA
Environmental	Temperature	Operation: +5 – +40 °C Storage: -25 – +70 °C Transport: -25 – +70 °C
	Humidity	Operation: 15 % – 93 % non-condensing Storage: < 93 % non-condensing
	Atmospheric	700 hPa – 1060 hPa
IP classification		IP42
Degree of safety in presence of inflammable anesthetics		The device is not intended for use with flammable anaesthetic gases
Applied part		Mattress

<b>MATTRESS SPECIFICATION</b>		
Product	Size	Weight
CuroCell® Ci10	80/85/90/100/105/120 x 200 x 10 cm	3,7 kg (80x200 cm)
CuroCell® Ci10 PRO	80/85/90/100/105/120 x 200/210 x 10 cm	5,2 kg (80x200 cm)
CuroCell® Ci17	80/85/90/100/105/120 x 200 x 17 cm	10,0 kg (80x200 cm)
CuroCell® Ci17 PRO	80/85/90/100/105/120 x 200/210 x 17 cm	11,5 kg (80x200 cm)
CuroCell® Ci20	80/85/90/100/105/120 x 20 cm	11,0 kg (80x200 cm)
CuroCell® Ci20 PRO	80/85/90/100/105/120 x 200/210 x 20 cm	12,5 kg (80x200 cm)




























**Note:** Care of Sweden reserves the right to modify the product specification at any time.









## 11.1 Standards

The system is tested and approved according to the following European standards where applicable requirements are met.

IEC 60601-1	EN ISO 10993	ISO 3746
IEC 60601-1-2	EN 12182	ISO 11201
IEC 60601-1-11	EN 597-1	
IEC 60601-1-6	EN 597-2	
IEC 62304	EN ISO 14971	

## 11.2 Symbol key

<b>Symbols to convey medical device information</b>			
	CE-marked in accordance with Medical Device Regulation (EU) 2017/745		Manufacturer
	Medical Device		Distributor
	UDI		Date of manufacture
<b>Symbols for traceability and product information</b>			
	Item number		LOT
	SN-number		Type BF
	IP class (Enclosure class)		Class II Equipment (double insulated). Indicated on the power supply.
<b>Symbols for patient information</b>			
XXXX-XX-XX	Year-Month-Day		Foot placement
	Pressure ulcer category		Recommended patient weight
	Counteracts shear		Do not rotate
	Heel function		Do not turn around
	Place on top of existing mattress		Placed directly on the bed base
	Do not place directly on bed		Do not place on top of another mattress
	Minimum length xxx cm		The mattress should be used with the patient lying lengthways
	Consult instructions for use		Refer to instruction manual/booklet

<b>Symbols for cleaning and recycling</b>			
	Machine wash at 95 °C		Drip dry
	Tumble dry		Do not iron
	Wipe clean		Do not dry clean
	Chlorine		Do not dispose of with household waste

## 12 Other information

### 12.1 Recommended life time

The estimated lifetime of this product is 5 years, provided that the Service and maintenance chart in section 9.1 is followed.

### 12.2 Disassembly and recycling

Except for certain parts of the control unit, energy recovery is possible for almost all material in CuroCell® products through incineration in waste incineration facilities.

#### Control unit:

The air tube connector (marked 'CPR') is easy to disassemble and is sorted as "plastic waste". The other parts of the control unit must not be disassembled and are sorted as "electronic waste".

#### Mattress:

A used CuroCell® mattress should be taken to a recycling center. The product is sorted as 'combustible waste'.

#### Note!

- If it is assessed that the product is or could be contaminated (e.g., used by patients with a known bloodborne infection), the product must be handled in accordance with the healthcare provider's or local authority's procedures for contaminated waste.

### 12.3 Returns

Contact Care of Sweden or your local distributor before returning the product.



Manufactured by

 care of sweden

SUPPORTING LIFE

**Contact:**

Phone: + 46 771 106 600

Fax: + 46 325 12840

Email: [export@careofsweden.se](mailto:export@careofsweden.se)

Internet: [www.careofsweden.com](http://www.careofsweden.com)

**Address:**

Care of Sweden AB

P.O. Box 146

SE-514 23 Tranemo

SWEDEN

**Visit address:**

Fabriksgatan 5A

SE-514 33 Tranemo

SWEDEN

**Cargo address:**

Byns väg 4A

SE-514 33 Tranemo

SWEDEN

**Distributed by:**