


# CuroCell® 4 AD

Comfortable alternating pressure mattress system



CuroCell® 4 AD with Top part Olivia/Bottom part CC

# CuroCell® 4 AD

The CuroCell® 4 AD is used as an aid in the prevention and treatment of pressure ulcers up to and including category 4<sup>(1)</sup>. It is a functional and comfortable mattress replacement system and is intended for all healthcare environments.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene, quality and safety<sup>(2)(4)</sup>.

All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(3)</sup>, which contributes to reduce the risk of skin maceration.

## REFERENCES

(1) European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019

(2) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.

(3) SS-EN ISO 15496:2004, DIN 53122-1.

(4) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.



Always read the instructions for use prior to usage of the product.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB. © Care of Sweden AB, 2021. All rights reserved. Last update: 2021-06-02.

FEATURES & CHARACTERISTICS	FACTS
<ul style="list-style-type: none"> <li>• CPR function – rapid deflation</li> <li>• Digital pump with built-in sensor</li> <li>• Easy handling</li> <li>• Handles for transport and moving</li> <li>• Replaceable air cells of polyurethane with Cell-in-Cell system</li> <li>• Silent and energy-saving pump</li> <li>• Soft, elastic, vapour permeable<sup>(3)</sup> hygienic cover</li> <li>• Sound and light notification for low pressure</li> <li>• Sound and light notification for power failure</li> <li>• Static function with automatic return</li> <li>• Transport function</li> <li>• Good cleaning properties</li> </ul>	<ul style="list-style-type: none"> <li>• Type: alternating pressure mattress</li> <li>• Height: 21 cm</li> <li>• Function: an aid for prevention and treatment up to and including category 4<sup>(1)</sup></li> <li>• Recommended user weight: up to 200 kg</li> <li>• Power supply: 230 V/50 Hz</li> <li>• Noise level of pump: 19 dBA</li> <li>• Warranty: 2 years</li> <li>• Cleaning of cover: wipe with cleaning agent and/or disinfectants.</li> <li>• Machine wash max 95 °C, tumble drying</li> </ul>
OPTIONAL	SIZES
<ul style="list-style-type: none"> <li>• Hygienic cover with welded seams</li> </ul>	80/85/90/105/120x200 cm and 80/85/90x210 cm. If you need any other sizes please contact your distributor.