

# Optimal Hygiene Pillow

Comfortable & durable pillows

We offer a wide range of accessories, including pillows and pillow covers that are designed to meet healthcare needs and are also suitable for home environments.

Optimal Hygiene Pillow is a pillow that consists of ball fibers that can easily be fluffed up to maintain the original shape of the pillow.

Optimal Sens is an anatomical sleeping pillow. It is made of visco elastic foam that conforms to the contours of the neck and head, providing good support and relief to the neck and head.

## Covers for comfort and hygiene

Optimal Hygiene Pillows comes with a liquidproof hygiene cover for easy cleaning.

The liquidproof hygiene cover is manufactured in a four-way stretch<sup>(1)(2)</sup> fabric to reduce the risk of shear forces, and is vapour-permeable<sup>(3)</sup> to lower the risk of skin maceration. The hygiene cover also features a liquidproof zipper.



CE



Technical specification	
Material	Optimal Hygiene Pillow & Pillow Case Polyester/cotton Optimal Sens & Pillow Case Viskoelastiskt skum, 50 kg/m <sup>3</sup> Pillow Case Polyurethane/polyester
Sizes	Optimal Hygiene Pillow & Pillow Case 35x45 cm 50x60 cm 50x70 cm Optimal Sens & Pillow Case 30x48 cm Optimal Pillow Case cover 35x45 cm 50x60 cm 50x70 cm
Cleaning instructions	Cleaning Optimal Hygiene Pillow Machine wash at max 60 °C Cleaning Optimal Sens Spot cleaning with mild detergent and water. Cleaning Optimal Pillow case Wipe clean with detergent and/or disinfectant. Machine wash max. 95°C
CE-marking	Registered and marked in accordance with MDR (EU) 2017/745

## References

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019.
- (2) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements.
- (3) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instructions for use prior to use.

Care of Sweden and Curera are trademarks of Care of Sweden AB. © Care of Sweden AB, 2024. All rights reserved. Last update: 2024-01-22.