

Mattress covers

Covers for comfort and hygiene

Our support surfaces are supplied with a removeable and liquidproof hygiene cover for easy cleaning. Our hygiene covers are manufactured in a four-way stretch⁽¹⁾⁽³⁾ fabric to reduce the risk of shear forces, and is vapour-permeable⁽²⁾ to lower the risk of skin maceration.

Hygiene covers (available as full cover and as top part)



Olivia

- Color: grey
- Sewn seams
- Drip-proof zipper on three sides
- Material: 55 % polyester, 45 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2



Olivia H

- Color: grey
- Sewn seams
- Drip-proof zipper on three sides
- Handles on one side for moving the mattress
- Material: 55 % polyester, 45 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2



Marin

- Color: blue
- Sewn seams
- Zipper on one or three sides
- Material: 37 % polyester, 63 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2



Marin FR

- Color: blue
- Sewn seams
- Zipper on one short side
- Material: 55% polyester, 45 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2, SS 876 00 10 (high risk)



Steel

- Color: blue
- Welded seams
- Drip-proof zipper on four sides
- Material: 44 % polyamide, 56 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2



Stone

- Color: grey
- Welded seams
- Drip-proof zipper on four sides
- Material: 61 % polyamide, 39 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2, BS 6807 Crib 5

Bottom parts (used in combination with a top part)



Bottom part Evac

- Color: black
- Zippers on four sides
- Handles on short and long sides for moving the mattress
- Velcro straps for securing patients in emergency situations
- Material: 100 % polyester polyurethane
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2, BS 6807 Crib 5



Bottom part CuroCell®

- Color: black
- Zippers on four sides
- Material: 100 % polyester polyurethane
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2, BS 6807 Crib 5



Bottom part Duro

- Color: blue
- Zippers on four sides
- Material: 44 % polyamide, 56 % polyurethane coating
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2



Bottom part Optimal 5zon® Plus

- Color: black
- Zippers on four sides
- Handles on short and long sides for moving the mattress
- Material: 100 % polyester polyurethane
- Fire standard: SS 876 00 01 (normal risk), EN 597-1, -2

References

- (1) CDC (2003), Guidelines for Environmental Infection Control in Health-Care facilities, updated 2019
- (2) SS-EN ISO 15496:2004, DIN 53.122.
- (3) SS 876 00 20:2017, Healthcare textiles - Mattresses - Specifications and requirements

Technical specification

Cleaning instruction	Cleaning of covers: wipe with cleaning agent and/or disinfectants. Machine wash max 95 °C, tumble drying
----------------------	---



Always read the instructions for use prior to use.

Care of Sweden and CuroCell are trademarks of Care of Sweden AB.
© Care of Sweden AB, 2024. All rights reserved. Last update: 2024-05-23.